

TrueLine

TrueLine™ Centrifuge Tubes

Designed for Perfection



For more specific information on claims, visit the Certificates page at www.truelinelab.in.

Warranty/Disclaimer: Unless otherwise specified, all products are for research use only.

Superior Performance

- ▶ Improved maximum RCF
 - 15 mL rated to 17,000 x g
 - 50 mL rated to 20,000 x g
- ▶ Temperature range -80°C to 121°C
- ▶ 5-year shelf-life

Ease of Use

- ▶ Leakproof design
- ▶ Easy-on/easy-off cap
- ▶ Clear accurate graduations
- ▶ Large white marking spot

Reliable Quality

- ▶ US FDA approved virgin polypropylene
- ▶ RNase-/DNase-free
- ▶ Heavy-metal-free color concentrate
- ▶ Non-cytotoxic and free of trace metals
- ▶ Non-pyrogenic
- ▶ Sterile

Ordering Information

Cat. No.	Description	RCF	Qty/Pk	Qty/Cs
TLCT-15ML-S	15 mL centrifuge tubes	17,000	25	500
TLCT-50ML-S	50 mL centrifuge tubes	20,000	25	500

Corning Technologies India Pvt. Ltd.

2nd Floor, Pioneer Square, CRPF Road, Sector - 62
Gurugram, Haryana - 122 005, INDIA

Tel: +91 124 - 460 4000 fax: +91 124 - 460 4099

www.truelinelab.in

For listing of trademarks, visit www.truelinelab.in. All other trademarks are the property of their respective owners.