

TrueLine

TrueLine™ Universal Fit Pipet Tips

Designed for Precision



For more specific information on claims, visit the Certificates page at www.truelinelab.in.

Warranty/Disclaimer: Unless otherwise specified, all products are for research use only.

Superior Performance

- ▶ US FDA approved virgin polypropylene
- ▶ RNase-/DNase-free
- ▶ Endotoxin-free
- ▶ Human DNA-free, PCR inhibitor-free
- ▶ Non-cytotoxic and free of trace metals
- ▶ Non-pyrogenic
- ▶ Flawless 'super smooth' internal surface
- ▶ Unique orifice design prevents DNA denaturalization
- ▶ Graduation marks for easy and error-free visual checks
- ▶ Universal fit pipet tips

Ordering Information

Cat. No.	Description	Quantity
TLT-1000-B	1,000 µL pipet tips, blue, universal fit, beveled, bulk	500/Pack, 10 Packs/Case
TLT-1000-C	1,000 µL pipet tips, clear, universal fit, beveled, bulk	500/Pack, 10 Packs/Case
TLT-200-Y	200 µL pipet tips, yellow, universal fit, beveled, bulk	1,000/Pack, 20 Packs/Case
TLT-200-C	200 µL pipet tips, clear, universal fit, beveled, bulk	1,000/Pack, 20 Packs/Case

Corning Technologies India Pvt. Ltd.

2nd Floor, Pioneer Square, CRPF Road, Sector - 62
Gurugram, Haryana - 122 005, INDIA

Tel: +91 124 - 460 4000 fax: +91 124 - 460 4099

www.truelinelab.in

For listing of trademarks, visit www.truelinelab.in. All other trademarks are the property of their respective owners.